




## SPECTROGREEN




# Analysis of 200 g/L NaCl-Solutions by ICP-OES with Single Radial and Dual-Side-On Interface Plasma Observation

## Introduction


The electrolysis of NaCl using direct current is one of the basic processes in the chemical industry. In this chlor-alkaline electrolysis chlorine, caustic soda and hydrogen are produced simultaneously:



Chlorine is produced at the anode:



At the cathode the reaction is:



Especially with modern membrane processes extremely low concentrations of Al, Ba, Ca, Mg and Sr at a  $\mu\text{g/kg}$ -level in the feed brine are required to protect the electrolysis cells. The limiting concentrations of impurities in

NaCl for a selection of elements are given in Table 1.

Table 1: Typical feed brine purity requirements for the electrolysis of sodium chloride (selection)

Element	Limiting Value ( $\mu\text{g/kg}$ )	Typical Concentration ( $\mu\text{g/kg}$ )
Ca + Mg	20	2-3
Al	50	10-20
Ba	100	
I	100	< 50
Ni	10	
Hg	100	
Br	50000	
Sr	400	
Fe	50	
SiO <sub>2</sub>	5000	
Sum: Cr, Mo, Cu, Co, Mn, Ni, Zn, Pb, As, Sb	200	

## Analysis of 200 g/L NaCl-Solutions by ICP-OES with Single Radial 2 and Dual-Side-On Interface Plasma Observation



The SPECTROGREEN with dual side-on and single radial plasma observation was investigated for the analysis of brines used for the electrolysis of alkaline chlorides. The trace elements were directly determined in 200 g/L NaCl solutions. With both techniques the required sensitivity for Ca, Mg, Sr, Ba, Al, Fe, Cu, Hg and P, needed to fulfill high purity specifications, can be achieved with one single method. Using the DSOI interface, the sensitivity could be increased by an average factor of 1.5 compared to single radial plasma observation.

This report describes the principle methodology for the analysis of brine solutions. Detection limits for the typical range of elements are presented.

### Experimental

#### Instrumentation

All measurements were performed with the SPECTROGREEN ICP optical emission spectrometer (SPECTRO Analytical Instruments, Kleve, Germany) with single radial and dual side-on plasma observation. The latter provides enhanced sensitivity compared to single radial plasma observation for many elements. In addition, it offers a high matrix compatibility, large linear dynamic range and excellent precision without the need to change the plasma observation mode during analysis. The SPECTROGREEN features a Paschen-Runge spectrometer mount, employing the proprietary Optimized Rowland Circle Alignment (ORCA polychromator) technique. Consisting of two hollow section cast shells, optimized small volume and 15 linear CMOS detectors, the wavelength range between 165 and 770 nm can be analyzed, allowing complete spectrum capture within 3 s. Due to the unique reprocessing capabilities of the system, a new measurement is not required even if additional elements or lines are to be determined at a later point in time.

The optic is hermetically sealed and filled with argon, continuously circulated through a filter, which absorbs oxygen, water vapor and other species. High optical transmission in the UV is achieved, allowing the determination of non-metals as well as the use of prominent and interference free lines in this region.

An air-cooled, 27.12 MHz, free running type LDMOS ICP-generator is installed, which ensures excellent stability of the forward power even in the case of rapidly changing sample loads. All relevant ICP operating parameters are software controlled, allowing easy selection of the optimum operating conditions.

For sample introduction, a Cross Flow nebulizer, a Scott type spray chamber and an argon humidifier were used. The ICP operating conditions are given in Table 2.

*Table 2: Typical ICP operating conditions for Dual-Side-On Interface and Single Radial Interface*

	Dual-Side-On Interface	Side-On Interface
<b>Power</b>	1150 W	1400 W
<b>Observation Mode</b>	Dual-Side-On	Single Radial
<b>Coolant flow</b>	14.0 L/min	
<b>Auxiliary flow</b>	2.0 L/min	
<b>Nebulizer flow</b>	0.90 L/min	0.85 L/min
<b>Plasma Torch</b>	Quartz, fixed, 3.0 mm Injector tube	
<b>Spray Chamber</b>	Double pass (Scott type)	
<b>Nebulizer</b>	Cross Flow	
<b>Sample aspiration rate</b>	2 mL/min	
<b>Replicate read time</b>	57 s per replicate	







### Calibration

The standards were prepared with a matrix concentration of 200 g/L NaCl of “suprapure” quality [1]. For calibration, a commercially available multi-element standard [2] was used. All standards were acidified with 1%  $\text{HNO}_3$  (v/v) [3] and 2 mg/L Sc [4] was added as an internal standard to all solutions. The concentrations of the resulting calibration standards are given in the following table 3.

Table 3: Calibration standards

Element	Std.1 [mg/L]	Std.2 [mg/L]	Std.3 [mg/L]	Std.4 [mg/L]
Al	0	0.1	0.5	2
As	0	0.1	0.5	2
B	0	0.1	0.5	2
Ba	0	0.1	0.5	2
Be	0	0.1	0.5	2
Ca	0	0.1	0.5	2
Cd	0	0.1	0.5	2
Co	0	0.1	0.5	2
Cr	0	0.1	0.5	2
Cu	0	0.1	0.5	2
Fe	0	0.1	0.5	2
Hg	0	0.1	0.5	2
K	0	0.5	2.5	10
Li	0	0.1	0.5	2

Element	Std.1 [mg/L]	Std.2 [mg/L]	Std.3 [mg/L]	Std.4 [mg/L]
Mg	0	0.1	0.5	2
Mn	0	0.1	0.5	2
Mo	0	0.1	0.5	2
Ni	0	0.1	0.5	2
P	0	0.5	2.5	10
Pb	0	0.1	0.5	2
Sc	2	2	2	2
Si	0	0.1	0.5	2
Sn	0	0.1	0.5	2
Sr	0	0.1	0.5	2
Tl	0	0.1	0.5	2
V	0	0.1	0.5	2
Zn	0	0.1	0.5	2

## Results

Table 4 shows the selected wavelengths and the limits of detection (LOD) achieved. The LODs were calculated according to the equation [5]:

$$\text{LOD} = 3 \text{ RSD}_b \cdot c / 100 \text{ SBR}$$

Where:

- RSD<sub>b</sub>: - relative standard deviation of  
10 replicates of the blank (in %)
- c: - concentration of the standard
- SBR: - signal to background ratio

Table 4: Typical Limits of Detection (LOD) for the selected lines with Single Radial and Dual-Side-On plasma observation.

Element	$\lambda$ [nm]	LOD (3 $\sigma$ ) [ $\mu\text{g/L}$ ] Dual-Side-On	LOD (3 $\sigma$ ) [ $\mu\text{g/L}$ ] Single Radial
Al	167.078	0.3	0.3
As	189.042	4.5	6.5
B	182.641	1.0	1.5
Ba	455.404	0.4	0.6
Be	313.042	0.2	0.3
Ca	393.336	0.15	0.2
Cd	214.438	0.4	0.45
Cd	226.502	0.7	0.9
Co	228.616	1.4	1.8
Cr	267.716	1.7	2.3
Cu	324.754	1.6	4.0
Fe	259.941	1.3	1.8
Hg	184.950	1.5	2.8
K	766.491	39	75

Element	$\lambda$ [nm]	LOD (3 $\sigma$ ) [ $\mu\text{g/L}$ ] Dual-Side-On	LOD (3 $\sigma$ ) [ $\mu\text{g/L}$ ] Single Radial
Li	670.780	1.0	3.5
Mg	279.553	0.08	0.1
Mn	257.611	0.2	0.3
Mo	202.095	1.5	1.8
Ni	231.604	2.4	3.0
P	177.495	4.0	7.0
P	178.287	6.0	8.0
Pb	220.353	9.0	13
Si	251.612	3.0	4.5
Sn	189.991	2.5	3.5
Sr	407.771	0.1	0.2
Tl	190.864	7.0	8.0
V	311.071	3.5	4.0
Zn	213.856	0.2	0.35



### Conclusions

The SPECTROGREEN with dual-side-on interface or single radial plasma observation offers a simple, fast, accurate, precise and cost-efficient method for the analysis of brine solutions. Using the DSOI interface, the sensitivity could be increased by an average factor 1.5 of compared to single radial plasma observation. In combination with an autosampler, the SPECTROGREEN can be fully automated. Independent from the number of lines and elements an analysis (including three replicates and pre-flush) can be performed in less than four minutes.

### References

- [1] NaCl Suprapur®, Merck, Darmstadt, Germany
- [2] Bernd Kraft GmbH, Duisburg, Germany
- [3] HNO<sub>3</sub> Suprapur®, 65%, Merck, Darmstadt, Germany
- [5] Inorganic Ventures, Christiansburg VA, USA
- [4] P. W. J. M. Boumans, Spectrochim. Acta 46B, 431 (1991)

**www.spectro.com**



#### GERMANY

SPECTRO Analytical Instruments GmbH  
Boschstrasse 10  
D-47533 Kleve  
Tel. +49.2821.892.0  
spectro.sales@ametek.com

#### U.S.A.

SPECTRO Analytical Instruments Inc.  
50 Fordham Rd  
Wilmington 01887, MA  
Tel. +1 800 548 5809  
+1 201 642 3000  
spectro-usa.sales@ametek.com

#### CHINA

AMETEK Commercial  
Enterprise (Shanghai) CO., LTD.  
Part A1, A4 2nd Floor Building No. 1 Plot Section  
No. 526 Fute 3rd Road East; Pilot Free Trade Zone  
200131 Shanghai  
Tel. +86.400.022.7699  
spectro-china.sales@ametek.com

#### Subsidiaries:

► **FRANCE:** Tel. +33.1.3068.8970, spectro-france.sales@ametek.com ► **GREAT BRITAIN:** Tel. +44.1162.462.950, spectro-uk.sales@ametek.com  
► **INDIA:** Tel. +91.22.6196.8200, sales.spectroindia@ametek.com ► **ITALY:** Tel. +39.02.94693.1, spectro-italy.sales@ametek.com  
► **JAPAN:** Tel. +81.3.6809.2405, spectro-japan.info@ametek.co.jp ► **SOUTH AFRICA:** Tel. +27.11.979.4241, spectro-za.sales@ametek.com

**SPECTRO** operates worldwide and is present in more than 50 countries. For SPECTRO near you, please visit [www.spectro.com/worldwide](http://www.spectro.com/worldwide)  
© 2019 AMETEK Inc., all rights reserved, subject to technical modifications • B-19, Rev. 3 • Photos: SPECTRO, Adobe Stock  
Registered trademarks of SPECTRO Analytical Instruments GmbH •  **SPECTRO**: USA (3,645,267); EU (005673694); "SPECTRO": EU (009693763);  
SPECTROGREEN: EU (017931732)